



Article Automated Radiosynthesis, Preliminary In Vitro/In Vivo Characterization of OncoFAP-Based Radiopharmaceuticals for Cancer Imaging and Therapy

Francesco Bartoli ¹^[b], Philip Elsinga ², Luiza Reali Nazario ², Aureliano Zana ³^[b], Andrea Galbiati ³^[b], Jacopo Millul ³, Francesca Migliorini ³, Samuele Cazzamalli ³, Dario Neri ^{4,5}, Riemer H. J. A. Slart ^{2,6}^[b] and Paola Anna Erba ^{1,2,*}^[b]

- ¹ Nuclear Medicine, Department of Translational Research and Advanced Technologies in Medicine and Surgery, University of Pisa and Azienda Ospedaliero Universitaria Pisana, 56126 Pisa, Italy; francesco.bartoli@med.unipi.it
- ² Medical Imaging Center, University Medical Center Groningen, University of Groningen, 9712 CP Groningen, The Netherlands; p.h.elsinga@umcg.nl (P.E.); l.reali.nazario@umcg.nl (L.R.N.); r.h.j.a.slart@umcg.nl (R.H.J.A.S.)
- ³ Philochem AG, R&D Department, Libernstrasse 3, CH-8112 Otelfingen, Switzerland; aureliano.zana@philochem.ch (A.Z.); andrea.galbiati@philochem.ch (A.G.); jacopo.millul@philochem.ch (J.M.); francesca.migliorini@philochem.ch (F.M.); samuele.cazzamalli@philochem.ch (S.C.)
- Department of Chemistry and Applied Biosciences, Swiss Federal Institute of Technology,
- CH-8093 Zurich, Switzerland; dario.neri@philogen.com
- Philogen S.p.A., 53100 Siena, Italy
- ⁶ Biomedical Photonic Imaging Group, Faculty of Science and Technology, University of Twente, 7522 NB Enschede, The Netherlands
- Correspondence: paola.erba@unipi.it

Abstract: FAP-targeted radiopharmaceuticals represent a breakthrough in cancer imaging and a viable option for therapeutic applications. OncoFAP is an ultra-high-affinity ligand of FAP with a dissociation constant of 680 pM. OncoFAP has been recently discovered and clinically validated for PET imaging procedures in patients with solid malignancies. While more and more clinical validation is becoming available, the need for scalable and robust procedures for the preparation of this new class of radiopharmaceuticals continues to increase. In this article, we present the development of automated radiolabeling procedures for the preparation of OncoFAP-based radiopharmaceuticals for cancer imaging and therapy. A new series of [68Ga]Ga-OncoFAP, [177Lu]Lu-OncoFAP and [18F]AIF-OncoFAP was produced with high radiochemical yields. Chemical and biochemical characterization after radiolabeling confirmed its excellent stability, retention of high affinity for FAP and absence of radiolysis by-products. The in vivo biodistribution of [18F]AIF-NOTA-OncoFAP, a candidate for PET imaging procedures in patients, was assessed in mice bearing FAP-positive solid tumors. The product showed rapid accumulation in solid tumors, with an average of 6.6% ID/g one hour after systemic administration and excellent tumor-to-healthy organs ratio. We have developed simple, quick, safe and robust synthetic procedures for the preparation of theranostic OncoFAP-compounds based on Gallium-68, Lutetium-177 and Fluorine-18 using the commercially available FASTlab synthesis module.

Keywords: ⁶⁸Ga-labeled; ¹⁷⁷Lu-labeled; Al[¹⁸F]F-labeled; OncoFAP; automated radiosynthesis; positron emission tomography; radiopharmaceuticals; tumors imaging; theranostics

1. Introduction

Small organic ligands that selectively bind to tumor-associated antigens are increasingly applied as targeting delivery vehicles of small bioactive payloads such as radionuclides [1,2]. Fibroblast activation protein (FAP, FAP- α) is a type-II transmembrane serine



Citation: Bartoli, F.; Elsinga, P.; Nazario, L.R.; Zana, A.; Galbiati, A.; Millul, J.; Migliorini, F.; Cazzamalli, S.; Neri, D.; Slart, R.H.J.A.; et al. Automated Radiosynthesis, Preliminary In Vitro/In Vivo Characterization of OncoFAP-Based Radiopharmaceuticals for Cancer Imaging and Therapy. *Pharmaceuticals* **2022**, *15*, 958. https://doi.org/10.3390/ ph15080958

Academic Editors: Aurélie Maisonial-Besset and Fabien Caillé

Received: 23 June 2022 Accepted: 29 July 2022 Published: 2 August 2022

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). protease overexpressed on the surface of cancer-associated fibroblasts (CAFs) in the tumor microenvironment (TME) in most epithelial malignancies [3]. FAP has recently emerged as a novel in vivo imaging biomarker with prognostic and theranostic potential. Favorable biodistribution profile and tumors-to-background ratios have been demonstrated for several FAP-based radiopharmaceuticals in development for PET/CT imaging. A wide range of cancer types has been imaged successfully with FAP-radioligands based on small molecules, particularly in those clinical conditions where imaging with [¹⁸F]FDG, the most extensively used PET radiopharmaceutical in oncology, has encountered limitations [4–8].

To the best of our knowledge, OncoFAP is the FAP small organic ligand with the highest-affinity reported to date, with a dissociation constant of 680 pM compared with that of FAPI-04 of 1020 pM [9]. OncoFAP-radio and OncoFAP-fluorophore conjugates exhibit rapid, selective and efficient tumor-targeting performance in murine models of cancer [9]. Extensive pre-clinical characterization of [⁶⁸Ga]Ga-DOTAGA-OncoFAP demonstrates its excellent imaging performance with similar targeting properties in a head-to-head comparison with [⁶⁸Ga]Ga-FAPI-46 [8]. First-in-human data have recently confirmed the pan-tumoral targeting potential of [⁶⁸Ga]Ga-DOTAGA-OncoFAP, with high uptake in primary tumors (breast cancer, colon cancer, fibrosarcoma, hepatocellular carcinoma), lymph node and distant metastases and a rapid clearance from healthy organs [8].

Given the high expression of FAP in cancer and the absence of significant uptake of FAP-targeting small molecule tracers in normal tissues, OncoFAP-based radioligands have inherent potential for theranostic applications. DOTAGA-OncoFAP can be radiolabeled with radioactive metals entailing therapeutic potential such as lutetium-177, the current most popular radionuclide for therapeutic applications for neuroendocrine tumors [10] and prostate cancer [11], leading to a new family of theranostic radiopharmaceuticals.

The interest of using DOTAGA for binding metal radionuclides is not limited to ⁶⁸Ga-labeled imaging agents. Indeed, DOTAGA can be used to complex several other interesting radiometals, such as lutetium-177, scandium-44, indium-111, zirconium-89 and actinium-225 [12–14]. In this work, we present the chemical, biological and pre-clinical characterization of a series of radiotheranostics OncoFAP derivatives, specifically DOTAGA-OncoFAP, NODAGA-OncoFAP and NOTA-OncoFAP with gallium-68 and fluorine-18 for PET/CT imaging and with lutetium-177 for therapeutic applications. Further, the production of a simple preparation kit was also developed for [⁶⁸Ga]Ga-DOTAGA-OncoFAP and [⁶⁸Ga]Ga-NODAGA-OncoFAP radiotracers.

2. Results

OncoFAP-COOH was used as common precursor for the synthesis of DOTAGA-OncoFAP, NODAGA-OncoFAP and NOTA-OncoFAP (Figure 1). Final compounds were afforded by following standard coupling and deprotection procedures. Detailed synthetic procedures, chemical structures and product characterization are reported in the Supplementary Material.



Figure 1. Chemical structures of OncoFAP and its conjugates with radiometal chelators. OncoFAP is a portable small organic ligand that targets Fibroblast Activation Protein in solid tumors with ultra-high affinity. Structures of OncoFAP and its DOTAGA, NODAGA and NOTA conjugates are presented. DOTAGA-OncoFAP conjugate was radiolabeled with gallium-68 and lutetium-177. NODAGA-OncoFAP conjugate was radiolabeled with gallium-68, lutetium-177 and fluorine-18. NOTA-OncoFAP- derivative was radiolabeled with gallium-68 and fluorine-18.

2.1. Radiosynthesis

2.1.1. [68Ga]Ga-OncoFAP-Derivatives

The production of [68 Ga]Ga-DOTAGA-OncoFAP, [68 Ga]Ga-NODAGA-OncoFAP and [68 Ga]Ga-NOTA-OncoFAP using a 5–10 min automatic synthesis procedure at a reaction temperature of 95 °C was highly reproducible, with high radiochemical purity (RCP) and high radiochemical yields (RCY, Tables 1–4). The molar activity was 20–30 GBq/µmol with excellent RCY, exceeding 80% for >15 µg precursors amount. The best RCY were obtained using 25 µg of OncoFAP-derivates, a reaction time of 5 min and a reaction pH of 4.2. No significant differences were observed increasing the amount of precursor to 20 µg and 25 µg or increasing the reaction time to 10 min. The validation data of our synthesis conditions are reported in Table 5.

Table 1. Different synthetic conditions used for [68Ga]Ga-NOTA-OncoFAP.

Radiosynthesis [⁶⁸ Ga]Ga-NOTA-OncoFAP			
	Condition 1	Condition 2	Condition 3
NOTA-OncoFAP (µg)	15	20	25
Reaction time (min)	10	5	5
Starting activity (MBq)	715 ± 115	780 ± 150	770 ± 135
Final activity (MBq)	580 ± 110	641 ± 135	650 ± 116
RCY (%)	78 ± 3	82 ± 2	84 ± 1
RCY decay corrected (%)	89 ± 3	88 ± 2	90 ± 1
Final molar activity (GBq/µmol)	39 ± 7	32 ± 7	26 ± 5
RCP (%)	99.1 ± 0.1	99.1 ± 0.1	99.1 ± 0.1

Radiosynthesis [⁶⁸ Ga]Ga-NODAGA-OncoFAP				
	Condition 1	Condition 2	Condition 3	
NODAGA-OncoFAP (µg)	15	20	25	
Reaction time (min)	10	5	5	
Starting activity (MBq)	1085 ± 95	1000 ± 74	1100 ± 60	
Final activity (MBq)	880 ± 70	840 ± 56	954 ± 28	
RCY (%)	81 ± 2	84 ± 2	87 ± 2	
RCY decay corrected (%)	93 ± 2	90 ± 2	93 ± 3	
Final molar activity (GBq/µmol)	59 ± 5	42 ± 3	38 ± 2	
RCP (%)	99.2 ± 0.1	99.2 ± 0.1	99.2 ± 0.1	

 Table 2. Different synthetic conditions used for [⁶⁸Ga]Ga-NODAGA-OncoFAP.

 Table 3. Different synthetic conditions used for [68Ga]Ga-DOTAGA-OncoFAP.

Radiosynthesis [⁶⁸ Ga]Ga-DOTAGA-OncoFAP			
	Condition 1	Condition 2	Condition 3
DOTAGA-OncoFAP (µg)	15	20	25
Reaction time (min)	10	5	5
Starting activity (MBq)	1159 ± 93	1208 ± 93	1159 ± 57
Final activity (MBq)	924 ± 99	942 ± 60	915 ± 33
RCY (%)	75 ± 3	78 ± 2	79 ± 2
RCY decay corrected (%)	84.0 ± 3	87 ± 2	88 ± 2
Final molar activity (GBq/µmol)	55 ± 6	42 ± 3	33 ± 1
RCP (%)	99.3 ± 0.1	99.3 ± 0.1	99.4 ± 0.1

Table 4. Acceptance criteria and average result of five [⁶⁸Ga]-derivatives synthesis obtained with condition of condition 3.

Parameters	Acceptance Criteria	Average NOTA-OncoFAP	Average NODAGA-OncoFAP	Average DOTAGA-OncoFAP
[68Ga]-labeled activity MBq	200–1300	520–580	900–970	850-920
Volume (mL)	12	12	12	12
Aspect	Clear colorless solution	Clear colorless solution	Clear colorless solution	Clear colorless solution
Half-life (min)	62–74	66.8	66.9	67.1
- ⁶⁸ Ga ³⁺ - ⁶⁸ Ga colloidal - RCP	≤2% ≤3% ≥95%	$\begin{array}{c} 0.1\pm 0.1\%\\ 0.8\pm 0.1\%\\ 99.1\pm 0.1\%\end{array}$	$\begin{array}{c} 0.1\pm 0.1\%\\ 0.7\pm 0.1\%\\ 99.2\pm 0.1\%\end{array}$	$\begin{array}{c} 0.1 \pm 0.1\% \\ 0.5 \pm 0.1\% \\ 99.4 \pm 0.1\% \end{array}$
Ethanol concentration (v/v)	<10%	$4.9\%\pm0.5\%$	$4.8\%\pm0.5\%$	$4.9\%\pm0.5\%$
⁶⁸ Ge breakthrough	<10 ⁻³ %	$< 10^{-4}\%$	<10 ⁻⁴ %	<10 ⁻⁴ %
рН	4.0-8.0	5.0 ± 0.3	5.0 ± 0.3	5.0 ± 0.3
Bacterial endotoxins	<175 EU/V	<6 EU/V	<6 EU/V	<6 EU/V
Filter integrity (psi)	\geq 50	>50	>50	>50

Radiosynthesis [⁶⁸ Ga]Ga-NODAGA-OncoFAP (Simple Preparation)			
	Formate Buffer	Acetate Buffer	
NODAGA-OncoFAP (µg)	40	40	
Reaction pH	3.0 ± 0.1	4.4 ± 0.1	
Starting activity (GBq)	1.1 ± 0.2	1.1 ± 0.2	
Final molar activity (GBq/µmol)	26 ± 5	26 ± 5	
- ⁶⁸ Ga ³⁺ content - Colloidal ⁶⁸ Ga content - RCP	$3.4 \pm 2.6\%$ $1.5 \pm 0.2\%$ $95.1 \pm 2.5\%$	$\begin{array}{c} 7.9 \pm 1.7\% \\ 4.6 \pm 1.7\% \\ 87.5 \pm 2.8\% \end{array}$	

Table 5. Average value obtained for the synthesis of [⁶⁸Ga]Ga-NODAGA-OncoFAP.

For the OncoFAP kit, the best result was obtained with 10 min radiolabeling at a reaction temperature of 95 °C with 40 μ g of precursor and pH of 3.2 (Tables S2 and S3). The reduction of the pH to 3 using formate buffer was necessary to reduce the colloidal gallium impurity below 1.5 \pm 0.2%. In fact, the use of the same condition of the automatic synthesis procedure (acetate buffer at pH 4.2) resulted in RCP of 88 \pm 2.4% with significant amount of both free gallium-68 and colloidal gallium-68 (Tables 5 and 6).

Table 6. Average value obtained for the synthesis of [⁶⁸Ga]Ga-DOTAGA-OncoFAP.

Radiosynthesis [⁶⁸ Ga]Ga-DOTAGA-OncoFAP (Simple Preparation)			
	Formate Buffer	Acetate Buffer	
DOTAGA-OncoFAP (µg)	40	40	
Reaction pH	3.0 ± 0.1	4.4 ± 0.1	
Starting activity (GBq)	1.1 ± 0.2	1.1 ± 0.2	
Final molar activity (GBq/µmol)	26 ± 5	26 ± 5	
- ⁶⁸ Ga ³⁺ content - Colloidal ⁶⁸ Ga content - RCP	$\begin{array}{c} 3.1 \pm 1.1\% \\ 2.9 \pm 0.3\% \\ 94.0 \pm 1.4\% \end{array}$	$5.6 \pm 2.7\% \ 6.4 \pm 1.8\% \ 88.0 \pm 2.4\%$	

2.1.2. [¹⁸F]AlF-NOTA/NODAGA-OncoFAP

The radiosynthesis of [¹⁸F]AIF-NOTA-OncoFAP via [¹⁸F]AIF-approach resulted in a RCY of 20.9% with a RCP of 89.6% (Table 7). The synthesis of [¹⁸F]AIF-NODAGA-OncoFAP resulted in very low RCY and low RCP (Table 8). The best result for the NOTA-precursor was obtained using 250 μ g of precursor with a reaction time of 25 min and a reaction temperature of 95 °C (Table 9).

Table 7. Different conditions used for the synthesis of [¹⁸F]AlF-NOTA-OncoFAP.

Radiosynthesis of [¹⁸ F]AlF-NOTA-OncoFAP			
	Condition 1	Condition 2	Condition 3
NOTA-OncoFAP (µg)	200	250	300
Initial activity (MBq)	3970 ± 203	9250 ± 147	16280 ± 976
Product activity (MBq)	678 ± 33	1750 ± 100	2967 ± 25
RCY (%)	16.9 ± 0.5	18.9 ± 0.3	18.2 ± 1.2
RCY decay corrected (%)	18.7 ± 0.5	20.9 ± 0.3	20.2 ± 1.2
Final molar activity (GBq/µmol)	3.05 ± 0.3	6.6 ± 0.8	9.3 ± 0.7
RCP (%)	92.0 ± 0.5	89.6 ± 0.8	90.5 ± 0.9

Radiosynthesis of [¹⁸ F]AlF-NODAGA-OncoFAP						
	Condition 1 Condition 2 Condition 3					
NODAGA-OncoFAP (µg)	200	250	300			
Initial activity (MBq)	3219 ± 150	4255 ± 185	6660 ± 220			
Product activity (MBq)	78 ± 15	96 ± 21	137 ± 15			
RCY (%)	2.4 ± 0.5	2.3 ± 0.5	2.1 ± 0.3			
RCY decay corrected (%)	2.6 ± 0.5	2.5 ± 0.5	2.3 ± 0.3			
Final molar activity (GBq/µmol)	0.36 ± 0.10	0.35 ± 0.10	0.42 ± 0.13			
RCP (%)	60.3 ± 5.2	68.4 ± 5.8	67.5 ± 4.8			

Table 8. Different conditions used for the synthesis of [¹⁸F]AlF-NODAGA-OncoFAP.

Table 9. Acceptance criteria and average result of five [¹⁸F]AlF-OncoFAP synthesis obtained with condition of batch 3.

Parameters	Acceptance Criteria	Average NOTA-OncoFAP	Average NODAGA-OncoFAP
[¹⁸ F]-labeled activity	200–4000 MBq	2806–3205 MBq	210–230 MBq
Aspect	Clear colorless solution	Clear colorless solution	Clear colorless solution
RCY (%)	95%	20.2 ± 1.5	2.0 ± 0.7
- ¹⁸ F ⁻ content - Other - RCP	$\leq 2\% \leq 3\% \geq 95\%$	$\begin{array}{c} 4.1 \pm 0.2\% \\ 5.9 \pm 0.4\% \\ 90.6 \pm 0.6\% \end{array}$	$\begin{array}{c} 12.4 \pm 2.5\% \\ 27.3 \pm 1.4\% \\ 60.3 \pm 1.3\% \end{array}$
Ethanol concentration (v/v)	<10%	$5.6\pm0.5\%$	$5.4\pm0.5\%$
pН	4.0-8.0	5.5 ± 0.3	5.5 ± 0.3
Bacterial endotoxins	<175 EU/V	<6 EU/V	<6 EU/V
Filter integrity (psi)	\geq 50	>50	>50

2.1.3. [177Lu]Lu-DOTAGA-OncoFAP

For the radiosynthesis of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP, we obtained high molar activities of about 100 GBq/µmol with high RCY and RCP > 99.8% (Tables 10 and 11) using a reaction time of 30 min and a reaction temperature of 95 °C at a pH 4. We validated the radiolabeling method repeating the synthesis 3 times at the highest molar activity, simulating a typical production of a batch for PRRT i.e., 9.2 GBq of Lutetium-177 (Table 11). The RCY of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP strictly depends on the reaction volume and on the amount of precursor used in synthesis. The preparation was easily purified through C18 cartridge to remove the excess free Lutetium-177. No loss of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP was observed in this process (Figure S10 in the Supplementary Material).

Radiosynthesis of [¹⁷⁷ Lu]Lu-DOTAGA-OncoFAP				
	Condition 1	Condition 2	Condition 3	Condition 4
DOTAGA-OncoFAP (µg)	13.3 ± 0.25	5.1 ± 0.15	8.2 ± 0.50	5.3 ± 0.15
Reaction volume (mL)	2.5	2.5	2.5	1.5
Initial activity (MBq)	592 ± 9	870 ± 6	1047 ± 19	573 ± 6
Product activity (MBq)	577 ± 8	307 ± 5	670 ± 25	533 ± 6
RCY %	97.4 ± 0.4	35.3 ± 1.9	54.8 ± 1.3	92.9 ± 0.75
Gentisic acid (mg)	0	20	20	20
Final molar activity (GBq/µmol)	43 ± 2	60 ± 5	83 ± 3	100 ± 3
RCP (%)	> 99.8	>99.8	> 99.8	> 99.8

Table 10. Different conditions used for the synthesis of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP.

Table 11. Acceptance criteria and average result of three synthesis of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP condition 4.

[¹⁷⁷ Lu]Lu-DOTAGA-OncoFAP				
Parameters	Acceptance Criteria	Average DOTAGA-OncoFAP		
[¹⁷⁷ Lu]-labeled activity (GBq)	5–9	7.4–7.7		
Aspect	Clear colorless solution	Clear colorless solution		
- ¹⁷⁷ LuCl ₃ - Others - [¹⁷⁷ Lu]Lu-DOTAGA-OncoFAP-	$\leq 2\% \leq 3\% \geq 95\%$	$1.9 \pm 0.5\%$ - 98.1 $\pm 0.5\%$		
Ethanol concentration (v/v)	<10%	$3.2\pm0.9\%$		
рН	4.0-8.0	5.5 ± 0.3		
RCY		$95\pm2.6\%$		
Bacterial endotoxins	<175 EU/V	<6 EU/V		
Final molar activity (GBq/µmol)		105 ± 5		
Filter integrity (psi)	\geq 50	>50		

2.2. Preliminary In Vitro/In Vivo Characterization of OncoFA-Derivatives

The co-elution experiments showed that all the OncoFAP derivatives described in this work retained the ability to form stable complexes with recombinant hFAP (Figure 2 and Figure S10).

All [⁶⁸Ga]Ga-OncoFAP and [¹⁸F]AlF-OncoFAP radio-conjugates were highly stable (>99% of intact compound) in 0.9% saline solution and in human plasma at 37 °C for 2 h. On the other hand, a stability test of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP radiolabeled batch showed the presence of radiolysis of about 10% every 24 h, which was significantly reduced to about 2% within 8 days by the addition of 20 mg of gentisic acid as a radical scavenger (Figure 3b) [15,16]. All the radiolabeled preparations of OncoFAP derivatives exhibited high hydrophilicity. Details are given in Supplementary Material.

Cell binding studies showed a similar trend of binding on FAP-positive SK-RC-52 cells for all the OncoFAP radio-conjugates. All the compounds displayed an intense binding at early time point (10 min), followed by a progressive decrease over time. As shown in Figure 4, a major difference in the absolute binding value was found among 20% for [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP, 2.3% for [⁶⁸Ga]Ga-DOTAGA-OncoFAP and 0.25% for [¹⁸F]AIF-NOTA-OncoFAP. Results of the nonspecific binding on SK-R-RC-52 wild-type (SK-RC-52.wt) showed a very low absolute binding value for all the compounds (Figure 4).



Figure 2. Co-elution experiments performed with [⁶⁸Ga]Ga-DOTAGA-OncoFAP conjugate. [⁶⁸Ga]Ga-DOTAGA-OncoFAP forms a stable complex with recombinant hFAP. [⁶⁸Ga]Ga-DOTAGA-OncoFAP (500 kBq) and hFAP (2 μ M, 100 μ L) eluted on a PD-10 column.



Figure 3. In vitro radiopharmaceuticals' stability in saline or human plasma at 37 °C. Small aliquots were analyzed by HPLC (**a**) 30-60-90-120-150-180 min for [⁶⁸Ga]Ga-DOTAGA-OncoFAP and (**b**) 1–3–6–12–24–48–72–96 h and 7 days for [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP.



Figure 4. Cell binding assays in SK-RC-52.hFAP and SK-R-RC-52 wild-type (SK-RC-52.wt) cells (density of 1.5×10^5 cells each tube). (a) [⁶⁸Ga]Ga-DOTAGA-OncoFAP and [¹⁸F]AlF-NOTA-OncoFAP (10 µL, 0.1 MBq) cells were incubated at 37 °C for 10, 30, 60 and 120 min. (b) [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP (10 µL, 0.1 MBq) cells were incubated at 37 °C for 10, 120, 240, 360 min and 24 and 48 h.

In vivo biodistribution of [¹⁸F]AlF-NOTA-OncoFAP (500 nmol/kg, ~4 MBq/Kg, RCP 89%) in athymic Balb/c AnNRj-Foxn1 mice bearing subcutaneous HT-1080.hFAP fibrosarcoma showed selective accumulation in FAP-positive tumors (6.6% ID/g, 1 h post intravenous injection), with excellent selectivity against healthy organs. Tumor-to-blood and tumor-to-kidney ratios of 6.5-to-1 and 4.3-to-1, respectively, were observed at the 1-h time point (Figure 5).



Figure 5. (**A**) Quantitative in vivo biodistribution and (**B**) tumor-to-organ ratios of [¹⁸F]AlF-NOTA-OncoFAP in nude mice bearing subcutaneous HT-1080.hFAP tumors. [¹⁸F]AlF-NOTA-OncoFAP was intravenously injected at 500 nmol/Kg dose (3.85 MBq/Kg), animals were sacrificed 1-h post-administration and distribution of the radiopharmaceutical was assessed by gamma-counter.

3. Discussion

In recent years, the development of tumor-targeting ligands resulted in the discovery of a new generation of diagnostic and therapeutic products with high uptake in cancer lesions and low accumulation in healthy organs. Lutathera[®] and [¹⁷⁷Lu]Lu-PSMA-617 provide examples of radioligand therapeutics with proven clinical efficacy and limited systemic toxicity [17,18]. Radioligand diagnostic and therapeutic products targeting FAP in the tumor microenvironment represent a promising class of compounds with pan-tumoral applicability and high selectivity for tumor lesions [19,20].

FAP-targeting strategies have gained growing relevance in nuclear medicine for the development of radiolabeling. Although anti-FAP antibodies have been known since the 1990s [21–24], the discovery of small organic FAP ligands in the last few years represented a revolution in the field of nuclear medicine. Various FAP-targeting agents coupled to different radionuclide chelators have been recently developed and characterized for their high affinity and selectivity towards FAP-positive tumors. These conjugates exhibit rapid accumulation in cancer lesions and low uptake in healthy organs both in preclinical murine models and cancer patients. Among FAP-specific small organic ligands, OncoFAP is the compound with the highest affinity described so far (Kd = 680 pM) [9]. Preliminary data with [68Ga]Ga-DOTAGA-OncoFAP have demonstrated the feasibility of 68Ga-radiolabeling and highly favorable targeting properties in both small animal and patients with cancer, validating [68Ga]Ga-OncoFAP as a new powerful alternative to clinically established PET tracers [8]. BiOncoFAP, a bivalent derivative of OncoFAP, has been recently described as a novel compound with high and prolonged tumor uptake in murine models of cancer [25]. With the aim to identify a clinical [¹⁸F]-OncoFAP candidate, we generated NODAGA-OncoFAP and NOTA-OncoFAP as two preclinical prototypes that could be labeled both with Fluorine-18 and with Gallium-68. NODAGA-OncoFAP rapidly complexes Gallium-68 with excellent product stability, thus opening promising avenues to the development of simple kit preparations

The translation of radiopharmaceutical agents from preclinical development to clinical applications requires the implementation of standard synthetic and radiolabeling procedures [26]. In order to facilitate clinical development of novel FAP-radiopharmaceuticals, we have developed highly reproducible automated methodologies for the efficient radio-labeling of OncoFAP-derivatives (DOTAGA-OncoFAP, NODAGA-OncoFAP and NOTA-OncoFAP) with Gallium-68, Fluorine-18 and Lutetium-177 radioisotopes. The procedures described in this article are easy to implement, safe and robust. We based these procedures on the FASTLab automated module, a radiolabeling platform already broadly applied to produce other radiopharmaceuticals at the commercial scale [27–29]. With the aim to enable production of [⁶⁸Ga]-DOTAGA/NODAGA-OncoFAP radiopharmaceuticals also in small radiopharmacy sites with limited financial resources, we have also generated a

radiolabeling single-vial cold kit. This approach, already proposed for other products (Illuccix[™] and NETSPOT[®]; SOMAKIT TOC) [30–32], renders the production of ⁶⁸Ga-based products as easy as the one of ^{99m}Tc-based radiopharmaceuticals. The "OncoFAP-kit", which already contains all necessary items (buffers, ligands and excipients) to be combined with the radioisotopes, allows efficient production of [⁶⁸Ga]Ga-NODAGA-OncoFAP and [⁶⁸Ga]Ga-DOTAGA-OncoFAP for PET imaging applications. The reasons for the development of such a strategy are related to the need to simplify and facilitate, while maintaining high standard of quality, efficiency and reproducibility the preparation of ⁶⁸Ga imaging agents. Indeed, the traditional approach of automatic synthesis of ⁶⁸Ga-based radiopharmaceuticals generates substantial investments, i.e., hot cells and synthesis modules, quality control equipment, and highly qualified personnel who may not be easily accessible to small-scale radiopharmacies, thus hampering the widespread use of ⁶⁸Ga-PET imaging agents and its equal accessibility for patients throughout the world. Consequently, interest has grown on the development of cold kits for PET tracers. The main advantages of a kit-based process are lower investments, ease of use, absence of purification steps (only a 0.22 µm sterilizing filter is required at the generator outlet to ensure the sterility of the final formulation), absence of EtOH and the lower final volume at the cost of higher amounts of precursor required, thus resulting in lower molar activity. Further, the lack of purification requires high-quality radioactive eluate.

Very few data are currently available on the direct comparisons of the operator's radiation exposure when producing ⁶⁸Ga-radiopharmaceuticals using on automatic synthesis modules and simple kits. Frinde et al., by measuring the dose at the extremities via proximity meters fixed to the first phalanx of each middle finger for a 1.85 GBq generator at the date of calibration, reported 70 μ Sv for the left hand and 132 μ Sv for the right hand for labeling with an automatic synthesis module but 179 μ Sv and 152 μ Sv for kit-based preparation of [68Ga]Ga-DOTANOC [33], The measurements made by Kleynhans et al. during the synthesis of [68Ga]Ga-HBED-CC-PSMA-11 demonstrated a whole-body exposure significantly lower for automated synthesis with 2.05 \pm 0.99 μ Sv versus 14.32 \pm 5.3 μ Sv for the kit-based preparation [34] as a result of the different radiolabeling environment in which the procedures were performed: a hot cell with 50 mm lead shield for the automated synthesis as compared to a class II biosecurity cabinet shielded with a 3 mm lead layer, with a 10 mm lead shielded desktop screen for the kit-based preparation. When the same 50 mm lead hot cell was used, the whole-body exposure dropped to $2 \pm 0.5 \,\mu\text{Sv}$ also for the kit-based preparation [35] with a dose at the extremities of about 1.5 ± 0.4 mSv. Therefore, reducing the exposure for simple kit labeling is possible, but dependent on the environment of the radiolabeling and on the experience of the operator.

As an alternative to Gallium-68 for diagnostic applications, we also developed two novel OncoFAP-conjugates bearing NODAGA and NOTA radiometal chelators, which are known to form stable complexes with Alluminum-Fluorine-18 [36–40]. We performed the radiochemical synthesis of both [¹⁸F]AIF-NODAGA-OncoFAP and [¹⁸F]AIF-NOTA-OncoFAP. The radiolabeling of NODAGA-OncoFAP with Alluminum-Fluorine-18 resulted in very low RCY and RCP, probably due to the ability of the NODAGA chelator to form a stable neutral complex with the Al³⁺ ion, thus partially preventing the incorporation of the fluorine-18 as also reported by Liu et al. [41]. On the contrary, NOTA-OncoFAP, which devoid a carboxylic acid function, was efficiently radiolabeled with Aluminium-Fluorine-18, with a good purity and a RCY of approximately 20%. DOTAGA-OncoFAP can be labeled both with Gallium-68 and Lutetium-177, and thus can represent a "theranostic pair". This approach has been already successfully implemented for the diagnosis and therapy of Neuro-Endocrine Tumors (Lutathera[®]/NETSPOT[®]) [10] and prostate cancers ([⁶⁸Ga]Ga-PSMA-11/[¹⁷⁷Lu]Lu-PSMA-617) [11]. Our data prove the feasibility of producing high-quality [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP with molar activities of 100 GBq/µmol, which is almost double the molar activities currently routinely achieved with [¹⁷⁷Lu]Lu-PSMA-617 [42]. Incorporating the high activity of Lutetium-177 represents a critical parameter to obtain a therapeutic radiopharmaceutical (i.e., in the range of 7.4–14.8 GBq) containing

relatively low amounts of precursor, thus avoiding saturation of targets at the tumor site [9]. Such results serve as the basis to provide a fast track for the clinical application of therapeutic OncoFAP derivatives such as [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP, retaining selective uptake in the FAP-positive tumor already 10 min after injection (32% injected dose [ID]/g), with a sustained uptake (i.e., higher than 20% ID/g) in the lesions over the first 6 h time window and the more recently developed [¹⁷⁷Lu]Lu-DOTAGA-BiOncoFAP, exhibiting a more stable and prolonged tumor uptake than [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP (~20% ID/g vs. ~4% ID/g, at 24 h p.i., respectively) favorable tumor-to-organ ratios with low kidney uptake as well potent anti-tumor efficacy when administered at therapeutic doses in tumor-bearing mice [25].

All OncoFAP derivatives based on Gallium-68 and Aluminum-fluorine-18 studied in this work show excellent stability without the need of a stabilizer if not the small percentage of ethanol found in the preparation coming from purification process. On the contrary [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP- shows an intrinsic instability when there are high radioactive concentrations, therefore necessitating the addition of gentisic acid as a radical scavenger normally formed due to the β^- decay of Lutetium-177. All the OncoFAP derivatives also retained high affinity towards hFAP as shown from the coelution experiment of [68Ga] Ga-DOTAGA-OncoFAP with hFAP and by the cell binding experiments on FAP-positive SK-RC-52 cells. It is worth noticing that the cell binding experiments demonstrated a large difference in binding values among the compounds, varying from 0.25% for [¹⁸F]AIF-NOTA-OncoFAP to 2.5% for [⁶⁸Ga]Ga-DOTAGA-OncoFAP and 22% for [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP. We attributed this large difference as a consequence of the difference in the molar activity in the different preparations [43]. The excellent binding proprieties of OncoFAP derivatives were confirmed by the in vivo biodistribution studies in tumor-bearing mice, which confirm the favorable tumor-targeting performance of [¹⁸F]AIF-NOTA-OncoFAP, comparable to that reported for [⁶⁸Ga]/[¹⁷⁷Lu]-DOTAGA-OncoFAP [8,25]. [¹⁸F]AIF-NOTA-OncoFAP selectively accumulates in solid lesion, with a high tumor-to-background ratio at early time points (i.e., 1 h after systemic administration, the biodistribution time normally used for PET images), thus supporting further development of the molecule as diagnostic PET radiopharmaceutical in clinical practice. For comparison, [68Ga]Ga-DOTAGA-OncoFAP biodistribution demonstrated beneficial tracer kinetics and high uptake in murine FAP-expressing tumor models with high tumorto-blood ratios of 8.6 \pm 5.1 at 1 h and 38.1 \pm 33.1 at 3 h p.i. [8] Application of [¹⁸F]AlF-NOTA-OncoFAP as a PET/CT imaging agent may result in significant logistical and clinical advantages as a consequence of the imaging characteristic of the isotope. Fluorine-18 is characterized by ~97% β^+ emission, 635 keV maximum positron energy, short β^+ trajectory with mean positron range of 0.27 mm in soft tissue. Moreover, it can be produced in large scale in cyclotrons [44,45] with the possibility of a subsequent delivery to satellite sites other than the production site, as the half-life of the isotope is long enough to allow this strategy ($t_{1/2} = 109.8$ min). Translation in the clinical setting will be of significant impact also in terms of daily organization of the PET/CT schedule. In particular, the use of fluorine-18-based FAPI-derivatives allow an easier patient preparation (no need of fasting and maintaining low glucose level respect the use of FDG) and offer a longer imaging time window of approximately 30–180 min after tracer injection as compared to gallium-68 FAPI derivatives. In addition, OncoFAP-derivatives can be exploited as PET/CT imaging agents for a broader window of malignancies, which are not efficiently detected by other marketed radiopharmaceuticals, such as [¹⁸F]FDG, [¹⁸F]FLT, [¹⁸F]F-MCH and [¹⁸F]F-PSMA-1007 [1,46]. Examples are represented by oesophageal [47], liver [48] and pancreatic cancers [4], brain primary tumors and metastases [49] or head and neck cancers [50]. However, the results from the synthesis using the Aluminum-Fluorine-18 approach still warrant further development and optimization before the method can be efficiently translated into clinical use. In fact, progress is needed to increase the labeling yield and molar activities while minimizing the precursor, thus avoiding a post-radiolabeling purification step to obtain a product with high molar activity.

12 of 18

Our results provide the reference for a robust and efficient radiolabeling method for OncoFAP-radioligand products that are easy to implement in clinical practice in small-, medium- and large-sized hospital radiopharmacy. The different radiosynthetic strategies presented in this article for each radioisotope can be chosen on the basis of the specific needs of the various clinical centers. Automated production of [⁶⁸Ga]Ga-DOTAGA-OncoFAP, a clinically validated PET tracer [8], and of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP will facilitate the implementation of this theranostic pair in clinical practice. Efficient production of [¹⁸F]AIF-NOTA-OncoFAP will represent the basis for the implementation of this novel radiotracer for large-scale PET imaging applications.

4. Materials and Methods

4.1. Chemical Synthesis

Detailed chemical procedures and compound characterization for the synthesis of NOTA-OncoFAP and NODAGA-OncoFAP are reported in the Supplementary Material. DOTAGA-OncoFAP has been produced as previously described [9].

4.2. Radiochemistry

4.2.1. [⁶⁸Ga]GaOncoFAP-Derivatives Synthesis

[⁶⁸Ga]Ga-DOTAGA/NODAGA/NOTA-OncoFAP (15–25 μ g) were synthesized in an FASTlab synthesis module (GE Healthcare); the configuration of the cassette is presented in the Figure S3 of the Supplementary Material. The synthesis was carried out using ⁶⁸Ga (t_{1/2} = 68 min, β^+ = 89%, and EC = 11%) automatically eluted with 0.1 M HCl (4.5 mL, TRASIS ALLinONE reagent kit) from a 1.85 GBq (50 mCi) ⁶⁸Ge/⁶⁸Ga radionuclide generator (Eckert & Ziegler 1850 MBq, GalliaPharm, Radiopharma GmbH) without pre-purification. About 3.5 mL of the [⁶⁸Ga]GaCl₃ eluate and the precursor, dissolved in sodium acetate (750 μ L, 0.7 M TRASIS ALLinONE reagent kit), were transferred into the reaction vessel. After stirring for 5 min at 95 °C under gentle nitrogen flow, the reaction mixture was loaded for purification onto the preactivated C18 cartridge (Waters, Milford, MA, USA), washed with 0.9% NaCl solution (5 mL, TRASIS ALLinONE reagent kit) and eluted with a mixture of 700 μ L of absolute ethanol (TRASIS ALLinONE reagent kit) and 800 μ L of water for injection (WFI, Fresenius kabi). Then, the product was diluted with 0.9% NaCl solution to obtain the final formulation obtaining a volume of 12 mL.

For the synthesis of [⁶⁸Ga]Ga-DOTAGA-OncoFAP and [⁶⁸Ga]Ga-NODAGA-OncoFAP, a simple preformulated radiopharmaceutical kit (OncoFAP-kit) was developed. This novel kit is 'ready-to-use' for radiolabeling with ⁶⁸Ga eluted by commercially available ⁶⁸Ge/⁶⁸Ga generators, as in the SomaKit TOC [32]. To this aim, we develop an approach consisting of generator elution with HCl (5 mL, 0.1 M TRASIS ALLinONE reagent kit) directly in the vial containing the precursor; at the end of the elution either sodium acetate (900 μ L, 0.7 M TRASIS ALLinONE reagent kit) to have a reaction pH of 4.5 or sodium formate (600 μ L, 1 M BioULTRA Formic acid solution 1.0 M in H₂O, Sigma Aldrich Taufkirchen, Germany and Sodium hydroxide solution BioUltra, for molecular biology, 10 M in H₂O) to have a reaction pH of 3 is added. Once the generator has been eluted into the vial containing the precursor and the appropriate buffer has been added, the vial is then placed in a thermoblock preheated to 98 °C and heated for 10 min. The synthesis was repeated in different conditions, varying the amount of precursor (10–40 μ g) and the reaction pH (3.2–4.5).

4.2.2. [¹⁸F]AlF-OncoFAP

Radiolabeling of NODAGA-OncoFAP and NOTA-OncoFAP (200–300 μ g) with ¹⁸F was performed via aluminum (Al³⁺) [¹⁸F]fluoride complex. In addition, for this synthesis, we used the FASTlab synthesis module. The configuration of the cassette is detailed in Figure S5 of the Supplementary Material.

Fluorine-18 was transferred to the module and trapped on a preactivated Sep-Pak light Accel plus QMA cartridge (Cl⁻ form: Waters, Milford, MA, USA). The cartridge was

washed with 6 mL of water (HPCE grade, Sigma Aldrich) and subsequently eluted with 500 μ L of the eluent solution, composed of 250 μ L of 0.9% NaCl solution (99.999% trace metals basis NaCl, Sigma Aldrich) in water for injection (WFI, Fresenius kabi) and 250 μ L absolute ethanol (TRASIS ALLinONE reagent kit), into a 5 mL reactor vial prefilled with 25 μ L of 2 mM aluminum chloride solution (AlCl₃, anhydrous, powder, 99.999% trace metals basis, Sigma-Aldrich) in sodium acetate buffer (0.1 M, pH 4.1). After stirring for 5 min at room temperature under gentle nitrogen flow to allow the formation of [¹⁸F]AlF, the solution of the precursor (600 μ L of 350 μ g/mL NOTA-OncoFAP or NODAGA-OncoFAP in sodium acetate 0.1 M pH 4.5) was added to the reactor vial, which was sealed and heated for 10 min at 95 °C. Next, the mixture was cooled to 40 °C. The reaction was diluted with 3.5 mL of 0.9% NaCl solution, loaded for purification onto the preactivated C18 cartridge (Waters, Milford, MA, USA), washed with 5 mL of 0.9% NaCl solution and eluted with 1.5 mL of absolute EtOH (TRASIS ALLinONE reagent kit). The product is diluted with 0.9% NaCl solution to obtain the final formulation obtaining a final volume of 12 mL. A method using semi-preparative HPLC purification is given in the supplementary information.

4.2.3. [177Lu]Lu-OncoFAP

 $[^{177}$ Lu]Lu-DOTAGA-OncoFAP was synthesized in an FASTlab synthesis module (GE Healthcare), at different molar activities ranging from 44 to 105 GBq/µmol. The configuration of the cassette is presented in Figure S7 of the Supplementary Material. Carrier-free Lutetium-177 (EndolucinBeta[®] 40 GBq/mL—pharmaceutical precursor, solution) with a concentration activity of 37 MBq/µL in a volume of 0.5 mL acetate buffer at pH 4.5 is transferred to the reactor vial. The reaction was scaled down to an amount of DOTAGA-OncoFAP ranging from 5.1 to 13.3 µg, which was dissolved in acetate buffer (1 mL, 0.1 M pH 4.5 TRASIS ALLinONE reagent kit) and aspirated into the reactor vessel. The solution was stirred for 30 min at 90 °C under gentle nitrogen flow, then loaded for purification onto the preactivated C18 cartridge, washed with 5 mL of 0.9% NaCl solution and eluted with eluent solution made of 700 µL of absolute ethanol (TRASIS ALLinONE reagent kit) and 800 µL of water for injection (WFI, Fresenius kabi). The product was diluted with 0.9% NaCl solution to obtain the final formulation. Further, the radiolabeling was also performed by adding 20 mg of gentisic acid [16] to assess its effects on the radiolytic stability as radical scavenger.

4.2.4. Reproducibility of the Method

To test the robustness of the synthesis method, we varied the reaction time and the amount of precursor used in the tests. We repeated each experiment at least three times to test the reproducibility.

4.3. Quality Controls

Radiochemical, chemical and radionuclide purity, pH, half-life, residual organic solvents, filter integrity and endotoxin content were assessed for all the radiolabeled preparations of OncoFAP derivatives.

Radiochemical and chemical purity were determined by HPLC (Berthold HERM LB500 radio detector with Jasco MD-2010 DAD at 250 nm) using a RP-18 column (Phenomenex Luna[®] column 150 \times 3 mm 3 µm 100 Å) with a flow of 0.6 mL/min at 40 °C. As mobile phase, water + TFA 0.1% (v/v, phase A) and acetonitrile + TFA 0.1% (v/v, phase B) were used with a gradient 0 to 0.5 min 90% eluent A, 0.5 to 10 min linear gradient elution from 90% to 0% eluent A, 10 to 14.5 min 0% eluent A, 14.50 to 15 min from 0% to 90% eluent A. The OncoFAP derivatives complex with the non-radioactive isotopes (Gallium-69, Fluorine-19 and Lutetium-175) were used as reference standards. iTLC on silica gel using citrate as a buffer eluent (0.1 M, pH 5) was used to quantify gallium-68 hydroxide with (Rf of 0–0.1), the labeled product (Rf between 0.3 and 0.5) and free gallium-68 (Rf between 0.6 and 0.8) both for the automated synthesis and the simple kit formulation.

For the determination of radionuclide purity, a sample of final product with a known activity and volume was analyzed (energetic spectrum from 0 to 1800 KeV) after one day in case of ⁶⁸Ga-compounds (⁶⁸Ga T_{decay} > 20 half-lives) and 1.5 days for [¹⁸F]AlF-OncoFAP derivatives (¹⁸F T_{decay} > 20 half-lives). To this aim, we selected only radioactive impurities with a half-life longer than 2 h. The pH of the final product was measured with pH indicator strips (pH range 2.0 to 9.0, MColorpHastTM, Merck, Milano, Italy). Final activity of ⁶⁸Ga/Al¹⁸F-OncoFAP was measured with an appropriate interval between the three measurements of 5 min. Residual organic solvents were measured using a gas chromatography system with a macrogol 20,000 column (30 m × 0.53 mm × 1 µm, Varian, Santa Clara, CA, U.S.) following the pharmacopoeia 2.2.28 (inlet 140 °C temperature, isocratic column temperature 50 °C, flame ionization detection temperature 250 °C during 5 min). Retention time of 2.9 min for ethanol and 4.75 min for *n*-propanol were used as reference standard.

Bubble point test of the Millex-GS vented filter was performed closing the vented portion of the filter and gradually pushing down the plunger inside a syringe filled with air to increase the pressure on the pressure gauge. The pressure when a continuous stream of air bubbles that appeared out of the 0.22 μ m membrane filter was measured as the product-wetted bubble point value.

Endotoxin determination was conducted with a Endosafe PTS Reader[®] rapid test (Charles River).

4.4. In Vitro Stability and Lipophilicity

The in vitro stability of the radiopharmaceuticals preparations in saline solution was evaluated on the bulk solution i.e., with the highest radioactive concentration. In addition, the in vitro stability of radiopharmaceutical preparations was tested by incubating 10 MBq of each radiolabeled compounds in 1 mL of human plasma at 37 °C. Small aliquots of bulk solution were analyzed by HPLC and iTLC and 10 μ L of plasma solution were analyzed by iTLC at 30, 60, 90, 120, 150 and 180 min for [⁶⁸Ga]Ga-DOTAGA/NODAGA-OncoFAP derivatives, 50, 100, 150 and 200-min for [¹⁸F]AlF-NOTA/NODAGA-OncoFAP or 1, 3, 6, 12, 24, and 2, 3, 4, 7 and 9 days for [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP, as show in Figure 3.

Lipophilicity was determined by partitioning between *n*-octanol and saline solution (LogP) or PBS (LogD_{7.4}) at room temperature using the conventional shake-flask method. An aliquot of the formulated solution containing ~500 kBq of radiopharmaceuticals was added to a tube containing 6 mL of *n*-octanol/saline solution or *n*-octanol/PBS (1:1 v/v, LogP or LogD_{7.4}, respectively). The tubes were shaken for 20 min, followed by centrifugation (5000 × *g* for 5 min) and separation of the phases. Aliquots of 1 mL of the organic and the aqueous phases were taken, and the activity was measured using an automated gamma counter. The LogP or LogD_{7.4} were calculated as the decadic logarithm of [activity (cpm/mL) in *n*-octanol]/[(activity (cpm/mL) in saline solution or PBS].

4.5. Characterization of OncoFAP Binding to hFAP

To assess the binding properties of radiolabeled OncoFAP-derivatives preparations to hFAP, [⁶⁸Ga]/[¹⁸F]AlF/[¹⁷⁷Lu]-OncoFAP derivatives (500 kBq) were incubated with hFAP (2 μ M, 100 μ L) and then loaded on PD-10 columns, pre-equilibrated with running buffer (50 mM Tris, 100 mM NaCl, and 1 mM ethylenediaminetetraacetic acid [EDTA], pH = 7.4) and eluted with running buffer. Fractions of the flow through (500 μ L) were collected, and the radioactivity intensity was measured immediately on a gamma counter (λ = 511 KeV for gallium-68 and fluorine-18 and λ = 208 KeV for lutetium-177) [9].

4.6. Cell Binding

Cell binding assays were performed in triplicate on SK-RC-52.hFAP and SK-R-RC-52.wt cells (density of 1.5×10^5 cells each tube), growth in RPMI-1640 medium (Gibco) containing 10% Fetal Bovine Serum (Gibco) and 1% Antibiotic-Antimycotic (Gibco) at 37 °C in 5% carbon dioxide (CO₂). After adding [⁶⁸Ga]Ga-DOTAGA-OncoFAP, [¹⁷⁷Lu]Lu-

DOTAGA-OncoFAP or [¹⁸F]AIF-NOTA-OncoFAP (10 μ L, 0.1 MBq) to the culture medium, the cells were incubated at 37 °C for 10, 30, 60 min for ⁶⁸Ga-derivatives and 10, 30, 60, 120, 240 min for ¹⁸F-derivatives or 10, 30, 60, 120, 240, 360 min and 24 and 48 h for the ¹⁷⁷Lu-derivative. As a control, nonspecific binding determination was carried out by incubating with the OncoFAP derivatives complex with the non-radioactive isotopes reference standard (100 μ M) with the same cell lines for 5 min before incubation with the radiolabeled compounds. After incubation, the cell pellets were washed twice with 1 mL of cold PBS (pH 7.4), centrifuged (400× *g*, 5 min) and the pellet measured by an automatic gamma counter (WIZARD2, PerkinElmer, Waltham, MA, USA).

4.7. Animal Studies

4.7.1. Implantation of Subcutaneous Tumors

Tumor cells were grown to 80% confluence in Dulbecco's Modified Eagle Medium (DMEM) supplemented with 10% Fetal Bovine Serum (Thermofisher, Waltham, MA, USA) and 1% antibiotic-antimycotic (Gibco) and detached with Trypsin-EDTA 0.05%. HT-1080.hFAP, cells (FAP positive cells) were resuspended in Hanks' Balanced Salt Solution medium. Aliquots of 5×10^6 cells (100 µL of suspension) were injected subcutaneously in the right flanks of female athymic Balb/c AnNRj-Foxn1 mice (6 to 8 weeks of age, Janvier-labs, Le Genest-Saint-Isle, France).

4.7.2. Biodistribution Studies in Tumor Bearing Mice

Female athymic Balb/c AnNRj-Foxn1 mice (6 to 8 weeks of age) were implanted in the right flank with HT-1080.hFAP tumors as described above. Tumors were allowed to grow to an average volume of 250 mm³. Mice were randomized (n = 5 per group) and injected intravenously with radiolabeled preparations of [¹⁸F]AIF-NOTA-OncoFAP (10 nmol/mice; 77 kBq; RCP:89%). Mice were euthanized 1 h after the injection by CO₂ asphyxiation. Organs were extracted, weighted, and radioactivity was measured with a Packard Cobra Gamma Counter. Values are expressed as percent ID/g \pm SD. We did not repeat the biodistribution experiments using [⁶⁸Ga]Ga-DOTAGA-OncoFAP and [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP since they are already reported [8,9].

5. Conclusions

Efficient production of pharmaceutical-grade [⁶⁸Ga]Ga-OncoFAP and [¹⁸F]AIF-OncoFAP can be achieved through automatic radiolabeling procedures described in this article. [⁶⁸Ga]Ga-OncoFAP can be obtained with a simple kit-based approach, thus allowing small radiopharmacies to offer this new PET radiotracer to cancer patients with limited financial burden. Furthermore, we have established methodologies to efficiently produce an OncoFAP-based radiotheranostic pair in which [⁶⁸Ga]Ga-DOTAGA-OncoFAP is used as a companion diagnostic of the targeted therapeutic agent [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP. The procedures described in this article are being implemented in the clinical protocols of theranostic trials with OncoFAP-derivatives to enable imaging and therapy of patients with solid tumors.

Supplementary Materials: The following supporting information can be downloaded at: https: //www.mdpi.com/article/10.3390/ph15080958/s1, chemical synthesis, Figure S1. radiolabeling reaction scheme of NODAGA-OncoFAP with ¹⁷⁷Lu; Figure S2. flow chart for radiosynthesis of [⁶⁸Ga]Ga-OncoFAP; Figure S3. Layout of the cassette for radiosynthesis of [⁶⁸Ga]Ga-OncoFAP derivatives; Figure S4. Flow chart for radiosynthesis of [¹⁸F]AlF-OncoFAP; Figure S5. Layout of the cassette for radiosynthesis of [¹⁸F]AlF-NOTA-OncoFAP; Figure S6. Flow chart for radiosynthesis of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP; Figure S7. Layout of the cassette for radiosynthesis of [⁶⁸Ga]Ga-DOTAGA-OncoFAP; Figure S8. Radio-HPLC trace performed after radiosynthesis of [⁶⁸Ga]Ga-DOTAGA-OncoFAP; Figure S9. Radio-HPLC trace performed after radiosynthesis of [¹⁷⁷Lu]Lu-DOTAGA-OncoFAP; Figure S10. Co-elution experiment of radiolabelled preparations; Table S1. Different conditions used for the synthesis [⁶⁸Ga]Ga-NODAGA-OncoFAP; Table S2. Different conditions used for the synthesis [⁶⁸Ga]Ga-NODAGA-OncoFAP; Table S3. Different conditions used

16 of 18

for the synthesis [⁶⁸Ga]Ga-NODAGA-OncoFAP; Table S4. Different conditions used for the synthesis [⁶⁸Ga]Ga-NODAGA-OncoFAP; Table S5. LogP and LogD_{7.4}.

Author Contributions: Conceptualization, all authors; methodology, F.B., L.R.N., A.Z., A.G., J.M. and P.A.E.; validation, F.B., L.R.N.; formal analysis, F.B.; investigation, F.B., A.Z., A.G. and J.M.; resources, P.A.E. and D.N.; data curation, F.B.; writing—original draft preparation, F.B.; writing—review and editing, all authors; visualization, all authors; project administration, P.A.E. and S.C.; funding acquisition, P.A.E. and D.N. All authors have read and agreed to the published version of the manuscript.

Funding: This research was partially funded by Philogen.

Institutional Review Board Statement: Animal experiments were conducted in accordance with Swiss animal welfare laws and regulations under the license number ZH006/2021 granted by the Veterinäramt des Kantons Zürich.

Informed Consent Statement: Not applicable.

Data Availability Statement: Data is contained within the article and Supplementary Material.

Acknowledgments: D.N. that is a co-founder and shareholder of Philogen (www.philogen.com), a Swiss-Italian Biotech company that operates in the field of ligand-based drug delivery. A.Z., A.G., J.M., F.M. and S.C. are employees of Philochem AG, the company that owns and has patented OncoFAP.

Conflicts of Interest: The authors declare the following competing financial interests related to the content of the present work; P.A.E., FB have contract research with Philogen; P.A.E. is Editorial Board member of the EJNMI, EJHI, European Radiology Experimental, Clinical and Translational Imaging, QJNM and is Board member of EANM (Education Chair); P.H. is Editor-in-Chief of EJNMI Radiopharmacy and Chemistry; R.S. is Associated Editor for EJNMII and EJHI; D.N. is a cofounder and shareholder of Philogen (http://www.philogen.com/en/ accessed on 1 August 2022), a Swiss-Italian Biotech company that operates in the field of ligand-based pharmacodelivery; A.Z., A.G., J.M., F.M. and S.C. are employees of Philochem A.G., the daughter company of Philogen that owns and has patented OncoFAP and its derivatives.

References

- 1. Kratochwil, C.; Flechsig, P.; Lindner, T.; Abderrahim, L.; Altmann, A.; Mier, W.; Adeberg, S.; Rathke, H.; Röhrich, M.; Winter, H.; et al. Ga-FAPI PET/CT: Tracer Uptake in 28 Different Kinds of Cancer. J. Nucl. Med. 2019, 60, 801–805. [CrossRef]
- Kulterer, O.C.; Pfaff, S.; Wadsak, W.; Garstka, N.; Remzi, M.; Vraka, C.; Nics, L.; Mitterhauser, M.; Bootz, F.; Cazzamalli, S.; et al. A Microdosing Study with ^{99m}Tc-PHC-102 for the SPECT/CT Imaging of Primary and Metastatic Lesions in Renal Cell Carcinoma Patients. J. Nucl. Med. 2021, 62, 360–365. [CrossRef] [PubMed]
- 3. Kalluri, R. The biology and function of fibroblasts in cancer. Nat. Rev. Cancer 2016, 16, 582–598. [CrossRef] [PubMed]
- Pang, Y.; Zhao, L.; Shang, Q.; Meng, T.; Feng, L.; Wang, S.; Guo, P.; Wu, X.; Lin, Q.; Wu, H.; et al. Positron emission tomography and computed tomography with [⁶⁸Ga]Ga-fibroblast activation protein inhibitors improves tumor detection and staging in patients with pancreatic cancer. *Eur. J. Nucl. Med. Mol. Imaging* 2022, 49, 1322–1337. [CrossRef] [PubMed]
- Guo, W.; Pang, Y.; Yao, L.; Zhao, L.; Fan, C.; Ke, J.; Guo, P.; Hao, B.; Fu, H.; Xie, C.; et al. Imaging fibroblast activation protein in liver cancer: A single-center post hoc retrospective analysis to compare [⁶⁸Ga]Ga-FAPI-04 PET/CT versus MRI and [¹⁸F]-FDG PET/CT. *Eur. J. Nucl. Med. Mol. Imaging* 2021, *48*, 1604–1617. [CrossRef] [PubMed]
- Giesel, F.L.; Kratochwil, C.; Schlittenhardt, J.; Dendl, K.; Eiber, M.; Staudinger, F.; Kessler, L.; Fendler, W.P.; Lindner, T.; Koerber, S.A.; et al. Head-to-head intra-individual comparison of biodistribution and tumor uptake of ⁶⁸Ga-FAPI and ¹⁸F-FDG PET/CT in cancer patients. *Eur. J. Nucl. Med. Mol. Imaging* 2021, 48, 4377–4385. [CrossRef]
- Gündoğan, C.; Kömek, H.; Can, C.; Yildirim, Ö.; Kaplan, İ.; Erdur, E.; Poyraz, K.; Güzel, Y.; Oruç, Z.; Çakabay, B. Comparison of ¹⁸F-FDG PET/CT and ⁶⁸Ga-FAPI-04 PET/CT in the staging and restaging of gastric adenocarcinoma. *Nucl. Med. Commun.* 2022, 43, 64–72. [CrossRef]
- Backhaus, P.; Gierse, F.; Burg, M.C.; Büther, F.; Asmus, I.; Dorten, P.; Cufe, J.; Roll, W.; Neri, D.; Cazzamalli, S.; et al. Translational imaging of the fibroblast activation protein (FAP) using the new ligand [⁶⁸Ga]Ga-OncoFAP-DOTAGA. *Eur. J. Nucl. Med. Mol. Imaging* 2022, 49, 1822–1832. [CrossRef]
- Millul, J.; Bassi, G.; Mock, J.; Elsayed, A.; Pellegrino, C.; Zana, A.; Dakhel Plaza, S.; Nadal, L.; Gloger, A.; Schmidt, E.; et al. An ultra-high-affinity small organic ligand of fibroblast activation protein for tumor-targeting applications. *Proc. Natl. Acad. Sci. USA* 2021, 118, e2101852118. [CrossRef]
- Strosberg, J.; El-Haddad, G.; Wolin, E.; Hendifar, A.; Yao, J.; Chasen, B.; Mittra, E.; Kunz, P.L.; Kulke, M.H.; Jacene, H.; et al. Phase 3 Trial of [¹⁷⁷Lu]-Dotatate for Midgut Neuroendocrine Tumors. N. Engl. J. Med. 2017, 376, 125–135. [CrossRef]

- Sartor, O.; de Bono, J.; Chi, K.N.; Fizazi, K.; Herrmann, K.; Rahbar, K.; Tagawa, S.T.; Nordquist, L.T.; Vaishampayan, N.; El-Haddad, G.; et al. Lutetium-177-PSMA-617 for Metastatic Castration-Resistant Prostate Cancer. N. Engl. J. Med. 2021, 385, 1091–1103. [CrossRef] [PubMed]
- 12. Price, E.W.; Orvig, C. Matching chelators to radiometals for radiopharmaceuticals. *Chem. Soc. Rev.* 2014, 43, 260–290. [CrossRef] [PubMed]
- Kostelnik, T.I.; Orvig, C. Radioactive Main Group and Rare Earth Metals for Imaging and Therapy. *Chem. Rev.* 2019, 119, 902–956. [CrossRef]
- 14. Kubícek, V.; Havlícková, J.; Kotek, J.; Tircsó, G.; Hermann, P.; Tóth, E.; Lukes, I. Gallium(III) complexes of DOTA and DOTAmonoamide: Kinetic and thermodynamic studies. *Inorg. Chem.* **2010**, *49*, 10960–10969. [CrossRef]
- 15. Trindade, V.; Balter, H. Oxidant and Antioxidant Effects of Gentisic Acid in a [¹⁷⁷Lu]-Labelled Methionine-Containing Minigastrin Analogue. *Curr. Radiopharm.* **2020**, *13*, 107–119. [CrossRef] [PubMed]
- 16. Joshi, R.; Gangabhagirathi, R.; Venu, S.; Adhikari, S.; Mukherjee, T. Antioxidant activity and free radical scavenging reactions of gentisic acid: In-vitro and pulse radiolysis studies. *Free Radic. Res.* **2012**, *46*, 11–20. [CrossRef] [PubMed]
- Baumgarten, J.; Groener, D.; Nguyen Ngoc, C.; Mader, N.; Chaurasia, M.; Davis, K.; Wichert, J.; Chun, F.K.H.; Tselis, N.; Happel, C.; et al. Safety and Efficacy of 177-Lutetium-PSMA-617 Radioligand Therapy Shortly after Failing 223-Radium-Dichloride. *Cancers* 2022, 14, 557. [CrossRef]
- Claringbold, P.G.; Turner, J.H. NeuroEndocrine Tumor Therapy with Lutetium-177-octreotate and Everolimus (NETTLE): A Phase I Study. *Cancer Biother. Radiopharm.* 2015, 30, 261–269. [CrossRef]
- 19. Gentles, A.J.; Newman, A.M.; Liu, C.L.; Bratman, S.V.; Feng, W.; Kim, D.; Nair, V.S.; Xu, Y.; Khuong, A.; Hoang, C.D.; et al. The prognostic landscape of genes and infiltrating immune cells across human cancers. *Nat. Med.* **2015**, *21*, 938–945. [CrossRef]
- 20. Valkenburg, K.C.; de Groot, A.E.; Pienta, K.J. Targeting the tumour stroma to improve cancer therapy. *Nat. Rev. Clin. Oncol.* 2018, 15, 366–381. [CrossRef]
- Schmidt, A.; Müller, D.; Mersmann, M.; Wüest, T.; Gerlach, E.; Garin-Chesa, P.; Rettig, W.J.; Pfizenmaier, K.; Moosmayer, D. Generation of human high-affinity antibodies specific for the fibroblast activation protein by guided selection. *Eur. J. Biochem.* 2001, 268, 1730–1738. [CrossRef] [PubMed]
- Welt, S.; Divgi, C.R.; Scott, A.M.; Garin-Chesa, P.; Finn, R.D.; Graham, M.; Carswell, E.A.; Cohen, A.; Larson, S.M.; Old, L.J. Antibody targeting in metastatic colon cancer: A phase I study of monoclonal antibody F19 against a cell-surface protein of reactive tumor stromal fibroblasts. *J. Clin. Oncol.* 1994, *12*, 1193–1203. [CrossRef] [PubMed]
- Tanswell, P.; Garin-Chesa, P.; Rettig, W.J.; Welt, S.; Divgi, C.R.; Casper, E.S.; Finn, R.D.; Larson, S.M.; Old, L.J.; Scott, A.M. Population pharmacokinetics of antifibroblast activation protein monoclonal antibody F19 in cancer patients. *Br. J. Clin. Pharmacol.* 2001, *51*, 177–180. [CrossRef] [PubMed]
- Mersmann, M.; Schmidt, A.; Rippmann, J.F.; Wüest, T.; Brocks, B.; Rettig, W.J.; Garin-Chesa, P.; Pfizenmaier, K.; Moosmayer, D. Human antibody derivatives against the fibroblast activation protein for tumor stroma targeting of carcinomas. *Int. J. Cancer* 2001, 92, 240–248. [CrossRef]
- 25. Galbiati, A.; Zana, A.; Bocci, M.; Millul, J.; Elsayed, A.; Mock, J.; Neri, D.; Cazzamalli, S. A novel dimeric FAP-targeting small molecule-radio conjugate with high and prolonged tumour uptake. *J. Nucl. Med.* **2022**, *63*. [CrossRef]
- Kolenc Peitl, P.; Rangger, C.; Garnuszek, P.; Mikolajczak, R.; Hubalewska-Dydejczyk, A.; Maina, T.; Erba, P.; Decristoforo, C. Clinical translation of theranostic radiopharmaceuticals: Current regulatory status and recent examples. *J. Label. Comp. Radiopharm.* 2019, 62, 673–683. [CrossRef] [PubMed]
- Khan, I.; Berg, T.C.; Brown, J.; Bhalla, R.; Wilson, A.; Black, A.; McRobbie, G.; Nairne, J.; Olsson, A.; Trigg, W. Development of an automated, GMP compliant FASTlab[™] radiosynthesis of [¹⁸F]GE-179 for the clinical study of activated NMDA receptors. *J. Label. Comp. Radiopharm.* 2020, 63, 183–195. [CrossRef]
- Edwards, R.; Greenwood, H.E.; McRobbie, G.; Khan, I.; Witney, T.H. Robust and Facile Automated Radiosynthesis of [¹⁸F]FSPG on the GE FASTlab. *Mol. Imaging Biol.* 2021, 23, 854–864. [CrossRef] [PubMed]
- Jiang, H.; Jain, M.K.; Cai, H. HPLC-free and cassette-based nucleophilic production of [¹⁸F]FDOPA for clinical use. Am. J. Nucl. Med. Mol. Imaging 2021, 11, 290–299.
- Illucix Information. Available online: http://illuccixhcp.com/wp-content/uploads/illuccix-prescribing-information.pdf (accessed on 2 February 2022).
- 31. Netspot Information. Available online: https://www.rxlist.com/netspot-drug.htm (accessed on 2 February 2022).
- 32. Somakit-TOC Information. Available online: https://www.ema.europa.eu/en/medicines/human/EPAR/somakit-toc (accessed on 2 February 2022).
- 33. Lepareur, N. Cold Kit Labeling: The Future of ⁶⁸Ga Radiopharmaceuticals? Front. Med. (Lausanne) 2022, 9, 812050. [CrossRef]
- 34. Kleynhans, J.; Rubow, S.; le Roux, J.; Marjanovic-Painter, B.; Zeevaart, J.R.; Ebenhan, T. Production of [⁶⁸Ga]Ga-PSMA: Comparing a manual kit-based method with a module-based automated synthesis approach. *J. Label. Comp. Radiopharm.* **2020**, *63*, 553–563. [CrossRef]
- Migliari, S.; Sammartano, A.; Scarlattei, M.; Baldari, G.; Janota, B.; Bonadonna, R.C.; Ruffini, L. Feasibility of a Scale-down Production of [⁶⁸Ga]Ga-NODAGA-Exendin-4 in a Hospital Based Radiopharmacy. *Curr. Radiopharm.* 2022, 15, 63–75. [CrossRef] [PubMed]

- Smith, G.E.; Sladen, H.L.; Biagini, S.C.; Blower, P.J. Inorganic approaches for radiolabelling biomolecules with fluorine-18 for imaging with positron emission tomography. *Dalton Trans.* 2011, 40, 6196–6205. [CrossRef] [PubMed]
- 37. Bruce Martin, R. Ternary complexes of Al3+ and F- with a third ligand. Coord. Chem. Rev. 1996, 149, 23–32. [CrossRef]
- 38. Antonny, B.; Chabre, M. Characterization of the aluminum and beryllium fluoride species which activate transducin. Analysis of the binding and dissociation kinetics. *J. Biol. Chem.* **1992**, *267*, *67*10–6718. [CrossRef]
- Li, L. The biochemistry and physiology of metallic fluoride: Action, mechanism, and implications. *Crit. Rev. Oral Biol. Med.* 2003, 14, 100–114. [CrossRef]
- McBride, W.J.; Sharkey, R.M.; Karacay, H.; D'Souza, C.A.; Rossi, E.A.; Laverman, P.; Chang, C.H.; Boerman, O.C.; Goldenberg, D.M. A novel method of ¹⁸F radiolabeling for PET. J. Nucl. Med. 2009, 50, 991–998. [CrossRef]
- 41. Liu, Y.; Hu, X.; Liu, H.; Bu, L.; Ma, X.; Cheng, K.; Li, J.; Tian, M.; Zhang, H.; Cheng, Z. A comparative study of radiolabeled bombesin analogs for the PET imaging of prostate cancer. *J. Nucl. Med.* **2013**, *54*, 2132–2138. [CrossRef]
- Chakraborty, S.; Chakravarty, R.; Shetty, P.; Vimalnath, K.V.; Sen, I.B.; Dash, A. Prospects of medium specific activity ¹⁷⁷Lu in targeted therapy of prostate cancer using ¹⁷⁷Lu-labeled PSMA inhibitor. *J. Label. Comp. Radiopharm.* 2016, 59, 364–371. [CrossRef]
- Luurtsema, G.; Pichler, V.; Bongarzone, S.; Seimbille, Y.; Elsinga, P.; Gee, A.; Vercouillie, J. EANM guideline for harmonisation on molar activity or specific activity of radiopharmaceuticals: Impact on safety and imaging quality. *EJNMMI Radiopharm. Chem.* 2021, *6*, 34. [CrossRef]
- 44. Le Bars, D. Fluorine-18 and medical imaging: Radiopharmaceuticals for positron emission tomography. J. Fluor. Chem. 2006, 127, 1488–1493. [CrossRef]
- 45. Sanchez-Crespo, A. Comparison of Gallium-68 and Fluorine-18 imaging characteristics in positron emission tomography. *Appl. Radiat. Isot.* **2013**, *76*, 55–62. [CrossRef] [PubMed]
- Treglia, G.; Muoio, B.; Roustaei, H.; Kiamanesh, Z.; Aryana, K.; Sadeghi, R. Head-to-Head Comparison of Fibroblast Activation Protein Inhibitors (FAPI) Radiotracers versus [¹⁸F]F-FDG in Oncology: A Systematic Review. *Int. J. Mol. Sci.* 2021, 22, 11192. [CrossRef] [PubMed]
- Zhao, L.; Chen, S.; Lin, L.; Sun, L.; Wu, H.; Lin, Q.; Chen, H. [⁶⁸Ga]Ga-DOTA-FAPI-04 improves tumor staging and monitors early response to chemoradiotherapy in a patient with esophageal cancer. *Eur. J. Nucl. Med. Mol. Imaging* 2020, 47, 3188–3189. [CrossRef]
- Dendl, K.; Finck, R.; Giesel, F.L.; Kratochwil, C.; Lindner, T.; Mier, W.; Cardinale, J.; Kesch, C.; Röhrich, M.; Rathke, H.; et al. FAP imaging in rare cancer entities-first clinical experience in a broad spectrum of malignancies. *Eur. J. Nucl. Med. Mol. Imaging* 2022, 49, 721–731. [CrossRef] [PubMed]
- Giesel, F.L.; Heussel, C.P.; Lindner, T.; Röhrich, M.; Rathke, H.; Kauczor, H.U.; Debus, J.; Haberkorn, U.; Kratochwil, C. FAPI-PET/CT improves staging in a lung cancer patient with cerebral metastasis. *Eur. J. Nucl. Med. Mol. Imaging* 2019, 46, 1754–1755. [CrossRef]
- Syed, M.; Flechsig, P.; Liermann, J.; Windisch, P.; Staudinger, F.; Akbaba, S.; Koerber, S.A.; Freudlsperger, C.; Plinkert, P.K.; Debus, J.; et al. Fibroblast activation protein inhibitor (FAPI) PET for diagnostics and advanced targeted radiotherapy in head and neck cancers. *Eur. J. Nucl. Med. Mol. Imaging* 2020, 47, 2836–2845. [CrossRef]